

Separation of protein and peptide monomers from their dimers, tetramers, and multimers presents a serious challenge to the separations scientist. Size-exclusion chromatography (SEC) and Tangential-Flow Filtration (TFF) (U.S. Pat. Nos. 5,256,294 and 5,490,937) have been used for separating monomers from aggregates but have limitations. SEC can separate monomers from multimers, and in some cases monomers from dimers. The main limitations of SEC are 1) limited load volumes (typically 5% of the bed volume) requiring large columns or multiple cycles, 2) and load protein concentration (low concentration feed stocks require pre-concentration or multiple cycles on the column. Higher protein concentrations can be more viscous, thereby reducing the efficiency of the separation). Historically TFF can separate protein multimers that are ten-fold larger than the monomer. U.S. Pat. No. 5,256,294.

U.S. Pat. Nos. 4,228,154 and 5,250,663 disclose separations of albumin from mixtures. U.S. Pat. No. 4,228,154 describes use of both cation-exchange and anion-exchange chromatography steps for the purification, without separation of monomer from multimers.

There is a need for separating monomers from dimers and multimers that is satisfactory, requires the use of only one ion-exchange step, and does not have the limitations of SEC or TFF.

Summary of the Invention

Accordingly, this invention provides a method for separating a polypeptide monomer from a mixture comprising dimers and/or multimers, wherein the method comprises applying the mixture to either a cation-exchange or an anion-exchange chromatography resin in a buffer, wherein if the resin is cation-exchange, the pH of the buffer is about 4-7, and wherein if the resin is anion-exchange, the pH of the buffer is about 6-9, and eluting the mixture at a gradient of about 0-1 M of an elution salt, wherein the monomer is separated from the dimers and/or multimers present in the mixture.

In this study it is demonstrated that ion-exchange chromatography—either anion or cation—is an effective means to separate protein or polypeptide monomers from their dimers and/or multimers.

Separations are performed using either step or linear gradient elution. Ion exchange has several advantages over the SEC and TFF methods described above. First, separation is independent of polypeptide concentration in the load and therefore no pre-concentration is required. Second, resins can be loaded to greater than 30 mg polypeptide/mL resin and still achieve excellent separations. Third, ion-

exchange resins are inexpensive and easy to use. Typical separations achieve enrichment of monomer to greater than 99.5% purity and yields in excess of 90%.

Brief Description of the Drawings

Figures 1A and 1B show separation of U266 IgE monomer from dimers and multimers on a RESOURCE Q™ anion-exchange column. The column was equilibrated in 25 mM Tris/pH 8, and eluted with a gradient from 0 to 0.5 M sodium chloride over 10 column volumes. Fig. 1A is full-scale; Fig. 1B is a close-up view to show the dimers and multimers.

Figures 2A1, 2A2, 2B, and 2C show separation of anti-IgE monoclonal antibody monomer from dimers and multimers. Figs 2A1 and 2A2 were run on a RESOURCE Q™ anion-exchange column. Fig. 2A1 is full-scale; Fig. 2A2 is a close-up view to show the dimers and multimers. Fig. 2B is a run on Q-SEPHAROSE FAST-FLOW™ resin. Fig. 2C is a plot of monomer and dimer/multimer observed in fractions, where the open dots are monomer and the solid dots are dimer. The monomer and dimer/multimer were determined using a SUPERDEX 200 HR™ 10/30 analytical size-exclusion column (Pharmacia Biotech). In all cases the columns were equilibrated in 25 mM Tris/pH 8. The gradient used in the Fig. 2A panels was 0 to 0.5 M sodium chloride over 40 column volumes. The gradient used for Fig. 2B (Q-SEPHAROSE FAST-FLOW™) was 0.05 to 0.2 M NaCl over 10 column volumes.

Figures 3A-C show separation of BSA monomer and dimer on a RESOURCE Q™ anion-exchange column at pH 8. The column was equilibrated in 25 mM Tris/pH 8, and eluted with a gradient from 0.125 to 0.275 M sodium chloride over 40 column volumes. Fig. 3A is purified monomer, Fig. 3B is purified dimer, and Fig. 3C is a commercial preparation of BSA (Bayer) that contains both monomer and dimer.

Figures 4A-C show separation of BSA monomer and dimer on a RESOURCE Q™ anion-exchange column at pH 6. The column was equilibrated in 20 mM sodium phosphate/pH 6, and eluted with a linear gradient from 0 to 0.5 M sodium chloride over 10 column volumes. Fig. 4A is purified monomer, Fig. 4B is purified dimer, and Fig. 4C is a commercial preparation of BSA (Bayer) that contains both monomer and dimer.

Figures 5A and 5B show separation of anti-IgE monoclonal antibody monomer from dimers and multimers on a RESOURCE S™ cation-exchange column at pH 6. The column was equilibrated in 20

mM sodium phosphate/pH 6, and eluted with a linear gradient from 0 to 0.05 M sodium chloride over 30 column volumes. Fig. 5A is the chromatogram from the separation, and Fig. 5B is a plot of monomer and dimer/multimer observed in fractions using the same method described in Figure 2, where the open dots are monomer and the solid dots are dimer.

Figures 6A and 6B show separation of BSA monomer and dimer on a RESOURCE S™ cation-exchange column at pH 4.3. The column was equilibrated in 20 mM sodium acetate/pH 4.3, then eluted with a gradient from 0 to 1 M sodium chloride over 20 column volumes. Fig. 6A is purified monomer, and Fig. 6B is purified dimer.

Detailed Description of the Preferred Embodiments

Definitions

As used herein, "polypeptide" refers generally to peptides and proteins having more than about ten amino acids. Preferably, the polypeptides are "exogenous," meaning that they are "heterologous," *i.e.*, foreign to the host cell being utilized, such as a human protein produced by *E. coli*. However, they may also be derived from a native source in which they are present naturally.

Examples of mammalian polypeptides include molecules such as, *e.g.*, renin, a growth hormone, including human growth hormone; bovine growth hormone; growth hormone releasing factor; parathyroid hormone; thyroid stimulating hormone; lipoproteins; 1-antitrypsin; insulin A-chain; insulin B-chain; proinsulin; thrombopoietin; follicle stimulating hormone; calcitonin; luteinizing hormone; glucagon; clotting factors such as factor VIIIc, factor IX, tissue factor, and von Willebrands factor; anti-clotting factors such as Protein C; atrial natriuretic factor; lung surfactant; a plasminogen activator, such as urokinase or human urine or tissue-type plasminogen activator (t-PA); bombesin; thrombin; hemopoietic growth factor; tumor necrosis factor-alpha and -beta; enkephalinase; a serum albumin such as human serum albumin; mullerian-inhibiting substance; relaxin A-chain; relaxin B-chain; prorelaxin; mouse gonadotropin-associated peptide; a microbial protein, such as beta-lactamase; DNase; inhibin; activin; vascular endothelial growth factor (VEGF); receptors for hormones or growth factors; integrin; protein A or D; rheumatoid factors; a neurotrophic factor such as brain-derived neurotrophic factor (BDNF), neurotrophin-3, -4, -5, or -6 (NT-3, NT-4, NT-5, or NT-6), or a nerve growth factor such as NGF; cardiostrophins (cardiac hypertrophy factor) such as cardiostrophin-1 (CT-1); platelet-derived growth

factor (PDGF); fibroblast growth factor such as aFGF and bFGF; epidermal growth factor (EGF); transforming growth factor (TGF) such as TGF-alpha and TGF-beta, including TGF- 1, TGF- 2, TGF- 3, TGF- 4, or TGF- 5; insulin-like growth factor-I and -II (IGF-I and IGF-II); des(1-3)-IGF-I (brain IGF-I), insulin-like growth factor binding proteins; CD proteins such as CD-3, CD-4, CD-8, and CD-19; erythropoietin; osteoinductive factors; immunotoxins; a bone morphogenetic protein (BMP); an interferon such as interferon-alpha, -beta, and -gamma; serum albumin, such as human serum albumin (HSA) or bovine serum albumin (BSA); colony stimulating factors (CSFs), *e.g.*, M-CSF, GM-CSF, and G-CSF; interleukins (ILs), *e.g.*, IL-1 to IL-10; anti-HER-2 antibody; superoxide dismutase; T-cell receptors; surface membrane proteins; decay accelerating factor; viral antigen such as, for example, a portion of the AIDS envelope; transport proteins; homing receptors; addressins; regulatory proteins; antibodies; and fragments of any of the above-listed polypeptides.

The preferred polypeptides of interest are mammalian polypeptides. Examples of such mammalian polypeptides include enzymes, hormones, cytokines, albumins, chemokines, immunotoxins, viral components, antibodies, neurotrophins, and antigens. Suitable such polypeptides encompass polypeptides such as HSA, BSA, anti-IgE, anti-CD20, anti-IgG, t-PA, gp120, anti-CD11a, anti-CD18, anti-VEGF, VEGF, TGF-beta, activin, inhibin, anti-HER-2, DNase, IGF-I, IGF-II, brain IGF-I, growth hormone, relaxin chains, growth hormone releasing factor, insulin chains or pro-insulin, NGF, NT-3, BDNF, and urokinase. Particularly preferred mammalian polypeptides include, *e.g.*, t-PA, gp120 (IIIb), anti-HER-2, anti-CD11a, anti-CD18, anti-VEGF, VEGF, BSA, HSA, anti-CD20, anti-IgE, anti-IgG, DNase, IGF-I, IGF-II, TGF-beta, IGFBP-3, IGFBP-2, IGFBP-1, growth hormone, NGF, NT-3, NT-4, NT-5, and NT-6. The polypeptide is more preferably an antibody or a serum albumin, more preferably, anti-IgE, anti-IgG, anti-Her-2, anti-CD11a, anti-CD18, anti-CD20, anti-VEGF, BSA, or HSA.

For purposes herein, the "mixture" contains monomers and either dimers or multimers or both dimers and multimers. Typically, the mixture is a biological fluid, which denotes any fluid derived from or containing cells, cell components, or cell products. Biological fluids include, but are not limited to, fermentation broth, cell culture supernatants, cell lysates, cleared cell lysates, cell extracts, tissue extracts, blood, plasma, serum, sputum, semen, mucus, milk, and fractions thereof. This definition includes cell-conditioned culture medium, which denotes a nutrient medium in which cells have been cultured and which contains cell products.

For purposes herein, "ion-exchange chromatography resin" refers to chromatography medium for anion- or cation-exchange separation.

As used herein, "elution salt" refers to an alkaline earth, alkali metal, or ammonium salt, i.e., a salt having a cation from the alkaline earth or alkali metal elements or an ammonium cation and having an inorganic or organic (hydrocarbon-based) anion. Examples of such salts include sodium chloride, ammonium chloride, sodium citrate, potassium citrate, potassium chloride, magnesium chloride, calcium chloride, sodium phosphate, calcium phosphate, ammonium phosphate, magnesium phosphate, potassium phosphate, sodium sulfate, ammonium sulfate, potassium sulfate, magnesium sulfate, calcium sulfate, etc. Preferred salts herein are chlorides or sulfates. The most preferred salt herein is sodium chloride.

As used herein, "multimers" refer to n-mers where n is 3-10, i.e., polymers that are not dimers but exclude aggregates. In contrast to multimers, aggregates have a value for n of greater than 10, and/or a molecular weight of greater than 2 million daltons, and/or are species contained in the excluded volume of analytical size-exclusion chromatography columns such as SUPEROSE 6™ (Pharmacia).

Modes for Carrying out the Invention

This invention relates to a method of separating monomers of polypeptides from their dimers or multimers or both. The method involves placing the mixture of monomers and dimers and/or multimers, from whatever source, in an equilibration buffer at a pH in the range of about 4 and 9 depending on whether the resin used for chromatographic separation is a cation- or anion-exchange resin. The resulting mixture is loaded onto either a cation-exchange or anion-exchange chromatography resin to capture all the n-mers (monomers, dimers, trimers, tetramers, etc.) present in the mixture. For ion-exchange column chromatography, ligands of general affinity can be used to achieve the desired selectivities and binding properties. The loading takes place in a buffer at a pH of about 6-9 if the resin is anion-exchange and about 4-7 if the resin is cation-exchange. The exact pH will depend, for example, on the isoelectric point of the polypeptide.

If the resin is a cation-exchange resin, prior to loading the mixture, the matrix can be equilibrated using several column volumes of a dilute, weak acid (e.g., four column volumes of 20 mM acetic acid, pH 3, or of 20 mM phosphoric acid, pH about 2.8). Following equilibration, the mixture is added and the column can be washed one to several times, prior to elution of the mixture, also using a weak acid solution such as a weak acetic acid or phosphoric acid solution. The buffer used for this purpose depends

on, e.g., the polypeptide and the anionic or cationic nature of the resin. For anion-exchange, preferably the buffer is TRIS or phosphate buffer; for cation-exchange, the buffer is preferably acetate or phosphate buffer.

Ion-exchange chromatography is typically carried out at a temperature of about 18-25°C, preferably about 20°C (room temperature). The preferred column loading is about 1 ml resin per 20-30 mg total polypeptide.

Following adsorption of the n-mer molecules to the ion exchanger, the mixture is eluted by contacting the resin with an elution salt having an appropriate ionic strength to displace the monomer from the matrix. An elution salt gradient is used of about 0 to 1 M. The gradient may be linear or stepwise. Preferably the gradient is from about 0 to 500 mM elution salt, more preferably 50 to 200 mM elution salt, and most preferably, 0 to 50 mM elution salt. Preferably the elution salt is a sodium salt, such as sodium chloride, although other elution salts and concentration gradients, known to those of skill in the art, also find use herein. The quantity of elution buffer can vary widely and will generally be in the range of about 2 to 40 column volumes, preferably 10 to 40 column volumes. Following elution, the eluate can be assayed for total monomeric concentration.

Suitable cation-exchange resins herein include a wide variety of materials known in the art, including those capable of binding polypeptides over a wide pH range. For example, carboxymethylated, sulfonated, agarose-based, or polymeric polystyrene/divinyl benzene cation-exchange matrices are particularly preferred. Other useful matrix materials include, but are not limited to, cellulose matrices, such as fibrous, microgranular, and beaded matrices; dextran, polyacrylate, polyvinyl, polystyrene, silica, and polyether matrices; and composites. These matrices include, for example, CM52 CELLULOSE™ (Whatman, Inc.); S-HYPERD™ and CM SPHERODEX™ (Secpracor); SP SEPHAROSE FF™, DEAE SEPHAROSE FF™, CM-SEPHAROSE™, and RESOURCE S™ (Amersham Pharmacia Biotech AB); and JT BAKER CSx™ (J.T. Baker, Inc.), as well as those containing the functional ligand R-SO₃⁻, preferably sulfopropyl resins, such as TOYOPEARL SP550C™ (Tosohaas) and FRACTOGEL EMD™ SO₃⁻-650 (m) (Merck). Other suitable materials for use in cation-exchange chromatography are within the knowledge of those skilled in the art.

Anion-exchange chromatography is carried out using media appropriate therefor, as are known in the art. Suitable media include, e.g., polymeric polystyrene/divinyl benzene resins and agarose-based

resins, as well as agarose beads, dextran beads, polystyrene beads, media that comprise an insoluble, particulate support derivatized with tertiary or quaternary amino groups., and supports derivatized with trimethylaminoethyl groups. Examples of suitable such media include DE92™ (diethylaminoethyl cellulose, Whatman); DEAE-CELLULOSE™ (Sigma), BAKERBOND ABX 40 mu™ (J.T. Baker, Inc.);
 5 DEAE resins such as FRACTOGEL EMD DEAE-650™ (EM Separations), FRACTOGEL EMD TMAE-650 (S)™ (EM Science, Gibbstown, NJ), TSK gel DEAE-SPW™ (Tosohaas), DEAE-SEPHAROSE CL-6B™ and chelating SEPHAROSE™ (Amersham Pharmacia Biotech AB), DEAE MERE SEP. 1000™ (Millipore), and DEAE SPHERODEX™ (Sepracor); RESOURCE Q™ and Q SEPHAROSE™ (QSFF) (Amersham Pharmacia Biotech AB); MACRO-PEP Q™ (Bio-Rad Laboratories, Hercules, CA); Q-
 10 HYPERD™ (BioSeptra, Inc., Marlborough, Mass); and the like. Other suitable anion-exchange chromatography materials, as well as the selection and use of these materials for the present application, are conventional in the art.

Purified fractions of monomer obtained from the ion-exchange chromatography may be further processed by subjecting them to any appropriate technique designed for downstream processing and
 15 purification. This will depend largely on the type of polypeptide and its intended use. Only one ion-exchange step is necessary to effect the desired separation of monomer from dimers and/or multimers in a mixture, although the invention does not exclude using more such steps if desired in the upstream or downstream processing of the polypeptide.

The invention will be more fully understood by reference to the following examples. They
 20 should not, however, be construed as limiting the scope of the invention. All literature and patent citations herein are incorporated by reference.

EXAMPLE I

This example shows the separation of anti-IgE monomers and bovine serum albumin monomers
 25 from dimers and multimers.

MATERIALS AND METHODS:

Resins:

Pharmacia Q-SEPHAROSE FAST FLOW™. 4mL to 235L bed volumes evaluated

Pharmacia RESOURCE S and RESOURCE Q™: 1mL prepacked columns
JT Baker CSx™, 0.46 X 5cm, 5μ particles

Proteins:

5 A. Humanized anti-IgE monoclonal antibodies (IgG₁) available from Genentech, Inc.: pI ~ 7.5, designated as E25 and E26. WO 93/04173 published March 4, 1993 describes humanized anti-IgE antibodies wherein a murine antibody directed against human IgE (MaE11) was used to provide the CDR regions that were substituted into an IgG1 immunoglobulin framework (rhuMaE25). See also Cacia et al., Biochemistry, 35: 1897-1903 (1996) for studies and further descriptions of E-25.

10 B. Monoclonal anti-IgE antibody prepared from the culture supernatants of an immortalized human myeloma cell line U266B1 (ATCC TIB 196) using affinity chromatography purification on an isolated anti-IgE antibody (Genentech MAE1). Specifically, five BALB/c female mice, age six weeks, were immunized in their foot pads with 10 μg of purified IgE in Ribi's adjuvant. Subsequent injections were done in the same manner at one and three weeks after the initial immunizations. Three days after
15 the final injection, the inguinal and popliteal lymph nodes were removed and pooled, and a single cell suspension was made by passing the tissue through steel gauze. The cells were fused at a 4:1 ratio with mouse myeloma P3X63-Ag8.653 (ATCC CRL 1580) in high glucose (DMEM) containing 50% w/v polyethylene glycol 4000. Alternatively, the immunizations were done in a similar manner except that 30 μg of IgE per injection were used and IgE fragment 315-347 (Kabat) was assayed as a prefusion
20 boost; or injections were given subcutaneously in two doses of 100 μg and a final booster of 50 μg, and spleen cells were used for the fusions.

The fused cells were then plated at a density of 2x10⁵ per well in 96-well tissue culture plates. After 24 hours HAT selective medium hypoxanthine/aminopterin/thymidine, Sigma, #H0262) was added. Of 1440 wells plated, 365 contained growing cells after HAT selection.

25 Fifteen days after the fusion, supernatants were tested for the presence of antibodies specific for human IgE using an enzyme-linked immunosorbent assay (ELISA). The ELISA was performed as follows, with all incubations done at room temperature. Test plates (Nunc Immunoplate) were coated for 2 hours with rat anti-mouse IgG (Boehringer Mannheim, #605-500) at 1 μg/ml in 50 mM sodium carbonate buffer, pH 9.6, then blocked with 0.5% bovine serum albumin in phosphate buffered saline

(PBS) for 30 minutes, then washed four times with PBS containing 0.05% TWEEN 20™ (PBST). Test supernatants were added and incubated two hours with shaking, then washed four times with PBST. Human IgE (purified from U266 cells as described above) was added at 0.5 µg/ml and incubated for one hour with shaking, then washed four times in PBST. Horseradish-peroxidase-conjugated goat anti-human IgE (Kirkegaard & Perry Labs, #14-10-04, 0.5 mg/ml) was added at a 1:2500 dilution and incubated for one hour, then washed four times with PBST. The plates were developed by adding 100 µl/well of a solution containing 10 mg of o-phenylenediamine dihydrochloride (Sigma, #P8287) and 10 µl of a 30% hydrogen peroxide solution in 25 ml phosphate citrate buffer, pH 5.0, and incubating for 15 minutes. The reaction was stopped by adding 100 µl/well of 2.5 M sulfuric acid. Data were obtained by reading the plates in an automated ELISA plate reader at an absorbance of 490 nm. For one antibody, 365 supernatants were tested and 100 were specific for human IgE. Similar frequencies of IgE specificity were obtained when screening for the other antibodies.

C. Bovine serum Albumin: pI 4.7 and 4.9 (Radola, Biochim. Biophys. Acta, 295: 412-428 (1973))

•Bayer Corp. P/N 81-024-2, "Bovine Albumin, Sulfhydryl Modified" (BSA Mix, blocked)

•ICN Biomedical Inc. P/N 810013, "Albumin Bovine" (BSA Mix, native)

•BSA monomer and dimer prepared in house from Bayer BSA (BSA Monomer and BSA Dimer, respectively)

Chromatography Systems:

Hewlett-Packard 1090™ HPLC

Pharmacia™ FPLC

Detection at 215 or 280 nm

Buffers: (see Table I for details)

Purified water

Tris•HCl

Sodium acetate

Sodium chloride

Sodium phosphate

Sodium citrate and citric acid

Sample Preparation:

Samples were diluted with the buffer used for equilibration (indicated in Table I below) to assure pH and conductivity matched starting column conditions. All samples were 0.2- μ m filtered prior to loading.

Chromatography:

Samples were introduced to the column using either an automatic or manual injector. All runs were performed at room temperature. Fractions were collected manually or with a PHARMACIA FRAC 100™ collector.

Chromatographic separation performance was evaluated by comparing elution profiles of BSA stock reagent and purified BSA monomer and dimer; the same was done for IgE and the monoclonal antibodies (MAb). Separation of IgE and MAb from their dimers and multimers was further evaluated by analyzing elution fraction using analytical size-exclusion chromatography. Plots of MAb MW forms vs. Fraction number were created. Recovery of IgE and MAb was determined spectrophotometrically by measuring absorbance at 280nm.

RESULTS:

The results are summarized in Table I below.

Table I

MONOMER-DIMER/MULTIMER SEPARATIONS

RESIN	PROTEIN	PH	EQUILIBRATION BUFFER	ELUTION	COMMENTS
Anion-Exchange					
QSFF™	MAb	8	Tris-HCl	linear gradient: 0 to 500 mM NaCl	Good separation
				linear gradient: 50 to 200 mM NaCl	Best separation
				step gradients to 200, 175, 150, 125 mM NaCl	Separation works
Resource	MAb	8	Tris-HCl		

	Q TM				
	Resource Q TM	U266 IgE	8	Tris-HCl	Removed aggregates and multimers
5	Resource Q TM	BSA Monomer	8	Tris-HCl	linear gradient: 0 to 1M NaCl, Good separation
		BSA Dimer		linear gradients: 150 to 550, 250 to 550mM NaCl	Excellent separation
		BSA Mix, native			
		BSA Mix, blocked		step gradients: 0.3/0.6, 0.38/0.6, 0.4M/0.6M NaCl	Some separation but fine control required
10	Resource Q TM	BSA Monomer	6	sodium citrate	linear gradient: Does not bind in citrate pH 6
		BSA Dimer			"
		BSA Mix, blocked			"
15	Resource Q TM	BSA Monomer	6	sodium phosphate	linear gradient: 0 to 0.5M NaCl in 10 CVs good separations
		BSA Dimer			
		BSA Mix, blocked			
Cation Exchange					
20	Resource S TM	BSA Monomer	6	sodium citrate	linear gradient: 0 to 0.5M NaCl in 10 CVs does not bind in citrate pH 6
		BSA Dimer			"
		MAb			"
	Resource S TM	MAb	6	sodium phosphate	linear gradient: 0 to 0.05M NaCl in 20 CVs equivalent to Q separation
					loaded to 16.5 mg/mL
25	Resource S TM	BSA Monomer	4.6	NaOAc buffer	linear gradient: 0 to 1M NaCl/40 CVs proteins somewhat resolved
		BSA Dimer			
		BSA Mix, blocked			
		BSA Monomer	4.3	NaOAc buffer	linear gradient: 0 to 1M NaCl/20 CVs better resolution than pH 4.6
		BSA Dimer			
		BSA Mix, blocked			

5	JT Baker CSx™	BSA Monomer	4.6	NaOAc buffer	linear gradient: 0 to 1M NaCl/12 CVs	proteins somewhat resolved
		BSA Dimer				
	JT Baker CSx™	BSA Monomer	4.3	NaOAc buffer	linear gradient: 0 to 1M NaCl/12 CVs	proteins somewhat resolved
		BSA Dimer				

Separations were evaluated using polymeric polystyrene/divinyl benzene resins (RESOURCE Q and S™), a silica-based resin (JT BAKER CSX™), and an agarose-based resin (Q-SEPHAROSE FAST FLOW™; QSFF). While separations were accomplished using any of these resins, separations worked especially well on Q-SEPHAROSE FAST FLOW™, RESOURCE Q™, and RESOURCE S™. The separation of BSA monomer and dimer from both suppliers looked very similar, suggesting the “Sulfhydryl Modified” material from Bayer did not alter the protein such that the species were easier to separate. It can be seen that phosphate buffer at pH 6 worked well, but no protein bound to the cation- or anion-exchange columns when 20 mM citrate buffer at pH 6 was used as equilibration buffer. Citrate buffer would be expected to work for both anion- and cation-exchange at a lower concentration, e.g., about 5 mM.

Recovery of monomeric IgE and MAbs to IgE on anion-exchange resins was typically greater than 90% at greater than 99.5% purity. Figures 1A and 1B show anion-exchange (RESOURCE™ Q) chromatograms in the separation of IgE monomers from dimers and multimers. Figures 2A1 and 2A2 show anion-exchange (RESOURCE™ Q) chromatograms in the separation of anti-IgE MAb monomers from dimers and multimers. Fig. 2B shows an anion-exchange (Q-SEPHAROSE FAST-FLOW™) chromatogram in the separation of anti-IgE MAb monomers from dimers and multimers. SEC (SUPERDEX 200 HR 10/30™) was used as an analytical method to determine the amount of monomer and multimer in samples from the ion-exchange separation, and Figure 2C shows the SEC analysis of fractions from Figure 2B. Separation of BSA monomer from dimer was readily achieved on anion-exchange resins at pH 8 and pH 6. See Figures 3A-C and 4A-C for chromatograms in the separation of BSA monomers from dimers and multimers by anion-exchange (RESOURCE™ Q) at pH 8 (Tris buffer) and at pH 6 (phosphate buffer), respectively.

Recovery and purity of MAb monomer from the cation-exchange resin was comparable to that of the anion-exchange resin. Figures 5A-B show cation-exchange (RESOURCE™ S) chromatograms in the

separation of anti-IgE MAb monomers from dimers and multimers at pH 6 (phosphate buffer).

Separations of BSA on cation-exchange resins could be performed at pH 4.6 and 4.3, 4.3 being somewhat better. Figures 6A-B show cation-exchange (RESOURCETM S) chromatograms in the separation of BSA monomers from dimers and multimers at pH 4.3 (acetate buffer).

- 5 In summary, mixtures of polypeptide mers were subjected to cation- or anion-exchange chromatography using a variety of resins and under a variety of pH and elution salt conditions, and successful separation was achieved. Based on results from four proteins with basic and acidic isoelectric points (two IgG₁ MAbs, IgE and serum albumin), the method demonstrates general applicability to separation of polypeptide monomers from their dimers and multimers.